



# Force Sensor

PS-2104

Built-in thumbscrew for mounting on a support rod



Holes for mounting to a cart or an accessory bracket

## Sensor Specifications

|                        |  |
|------------------------|--|
| Sensor Range:          | $\pm 50$ newtons (N)   |
| Accuracy:              | 1 %  |
| Resolution:            | 0.03 newtons (N)   |
| Max. Sample Rate:      | 1,000 sps  |
| Default Sample Rate:   | 10 sps   |
| Over-limit Protection: | Prevents damage from forces greater than 50 N  |
| ZERO Button:           | Tares the output to zero newtons before each use. Always tare with the sensor in the orientation used during the experiment. |

## Force Quick Start

The PS-2104 Force Sensor measures force in newtons.

## Additional Equipment Needed

- PASPORT Link Device (USB Link, **Xplorer**, etc.)
- EZscreen or DataStudio™ software (version 1.5 or later)
- Hook and rubber bumper attachments (included)
- Thumbscrew (included) used to mount on cart or Accessory Bracket (CI-6545)

## Equipment Setup

1. Connect the PASPORT Link Device to a USB port on your computer or USB hub.
2. Connect the sensor plug to a PASPORT Link Device.
3. The software launches when it detects a PASPORT sensor. From the PASPORTAL screen, select a point of entry:
  - an activity in the Workbook window,
  - EZscreen, or
  - DataStudio.



2



1

3



Click the Start  
Button to  
Record Data

Toggle  
Data  
Runs


Save  
Data

Open  
Data

Print  
Graph

Exit to  
DataStudio

Quit  
EZscreen



### EZscreen Specifications

|                              |   |
|------------------------------|---|
| <b>EZscreen Range:</b>       | ±50 newtons (N)   |
| <b>Recording Time:</b>       | up to 120 seconds   |
| <b>Scale-to-Fit:</b>         | Double-click the Graph to scale data                                  |
| <b>Information Tool:</b>     | Drag cursor over graph to display X,Y coordinate and slope at a point |
| <b>Export to DataStudio:</b> | Click <b>Exit to DataStudio</b> button                                |

## Force EZscreen

### EZscreen Activity-Acceleration Due to Gravity

1. Mount the Force Sensor on a horizontal support rod with hook pointing down.
2. Press the **Zero** button on the Force Sensor to tare it.
3. Obtain a variety of objects of known mass, such as the Hooked Mass Set (SE-8759). Hang a mass from the Force Sensor hook.
4. Click the **Start** button and record data for 10 seconds.
5. Click the **Stop** button.
6. Repeat Steps 2–5 for each mass.
7. Using your data and the formula:

$$F = mg$$

(where **F** equals the force exerted by the mass, **m** equals the mass, and **g** is the acceleration due to gravity), calculate an average value for the acceleration due to gravity.

8. Imagine that you performed this experiment on the Moon (where **g** is  $1.63 \text{ m/s}^2$ ). How would this experiment change? Explain.

Note: Masses between 200–1000g yield best results.